

Instrument: CS844

Carbon and Sulfur Determination in Low Carbon Ferroalloys

LECO Corporation; Saint Joseph, Michigan USA

Introduction

Ferroalloys are alloys of iron that contain a high level of one or more other primary elements. The principle ferroalloys consist of silicon, manganese, and chromium that are used as vehicles to introduce the alloying element into the molten metal when making steel or cast iron. For example, silicon is used to deoxidize steel and as an alloying element in cast iron. Manganese is used as an alloying element and mitigates the harmful effects of sulfur in cast iron and steel. Chromium increases corrosion resistance in stainless steels. Since carbon is the most important alloying constituent in steel and cast iron production, and sulfur is a harmful contaminant that negatively affects the mechanical properties of steel and cast iron, the determination of carbon and sulfur levels in the ferroalloy feed stock is a critical quality control parameter.

Sample Preparation

Samples should be a uniform, representative, powder, or granular material.

Accessories

528-018 or 528-018HP Crucibles (previously heated*); 502-173 LECOCEL II HP; 502-231 HP Iron Chip Accelerator; 763-266 LECOCEL; 501-078 Iron Powder; 501-636-HAZ V₂O₅ Accelerator; 773-579 Metal Scoop; 761-929 Tongs

*For optimal precision, ceramic crucibles are heated in a muffle or tube furnace (such as a LECO TF-4) at 1350 °C for a minimum of 20 minutes, or at 1000 °C for 40 minutes. The crucibles are removed from the furnace, allowed to cool for 1 to 2 minutes, and then are transferred to a desiccator for storage. Crucibles should be reheated if not used within four hours. After baking, crucibles should only be handled using clean tongs.

Reference Materials

LCRM[®], LRM[®], NIST, or other suitable reference materials, such as ferroalloy and steel reference materials.

Method Selection

Different methods can be used for the analysis of carbon and sulfur in ferroalloy materials on the CS844.

Method 1 utilizes LECOCEL II and iron chip accelerators to facilitate combustion, without the use of hazardous materials.

Method 2 utilizes iron powder, vanadium pentoxide, and LECOCEL as accelerators to facilitate combustion. This accelerator combination works well for ferroalloys and may improve sulfur recovery and precision. Even though the carbon blank intensity for this method is considered high, the blank is consistent enough to be properly subtracted from the analysis results. Vanadium pentoxide (V₂O₅) is considered a hazardous material, therefore, proper precautions should be taken.

Method Parameters

Analysis Parameters

Purge Time	10 s
Analysis Delay	20 s
Sample Cool Time	0 s
Furnace Mode	Constant
Furnace Power	100%

Element Parameters

	Carbon	Sulfur
Integration Delay	0 s	0 s
Starting Baseline	2 s	2 s
Use Comparator	Yes	Yes
Comparator Level	0.30%	0.30%
Minimum Integration Time	50 s	50 s
Maximum Integration Time	80 s	80 s
Significant Digits	5	5

Procedure / Method 1

- Prepare the instrument as outlined in the operator's instruction manual.
- Determine the instrument blank.
 - Login a minimum of three Blank replicates.
 - Add ~1.2 g of 502-173 LECOCEL II HP and ~0.8 g of 502-231 HP Iron Chip Accelerator to a previously heated crucible.
 - Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
 - Repeat steps 2b through 2c a minimum of three times.
 - Set the blank by following the procedure outlined in the operator's instruction manual.
- Calibrate/Drift Correct.
 - Login a minimum of three Standard replicates for each calibration/drift reference material to be used.
 - Weigh ~0.25 g of a calibration/drift reference material into the previously heated crucible.
 - Enter the mass and sample identification into the appropriate replicate fields.
 - Add ~1.2 g of 502-173 LECOCEL II HP and ~0.8 g of 502-231 HP Iron Chip Accelerator on top of the reference material.
 - Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
 - Repeat steps 3b through 3e a minimum of three times for each calibration/drift reference material utilized.
 - Calibrate/drift correct by following the procedure outlined in the operator's instruction manual.
 - Verify the calibration by analyzing ~0.25 g of another suitable reference material following steps 3b through 3e and confirm that the results are within an acceptable tolerance.
- Sample Analysis.
 - Login a sample with an appropriate number of replicates.
 - Weigh ~0.25 g of sample into a previously heated crucible.

- c. Enter the mass and sample identification into the appropriate replicate fields.
- d. Add ~1.2 g of 502-173 LECOCEL II HP and ~0.8 g of 502-231 HP Iron Chip Accelerator on top of the sample.
- e. Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
- f. Repeat steps 4a through 4e as necessary.

Typical Results / Method 1*

Sample	Mass (g)	% Carbon	% Sulfur
BCS 203/5	0.2499	0.0431	0.0168
Low Carbon	0.2514	0.0419	0.0165
Ferro-Chromium	0.2511	0.0433	0.0170
0.043% C	0.2517	0.0420	0.0161
	0.2506	0.0443	0.0165
	Avg =	0.0429	0.0166
	s =	0.0010	0.0004
Ferro-Chromium	0.2507	0.0325	0.0010
	0.2515	0.0323	0.0009
	0.2535	0.0322	0.0010
	0.2530	0.0332	0.0009
	0.2556	0.0340	0.0009
	Avg =	0.0329	0.0009
	s =	0.0008	< 0.0001
EURO 584-1	0.2541	0.0424	0.0296
Ferro-Titanium	0.2532	0.0419	0.0300
0.0445 ± 0.0043% C	0.2523	0.0430	0.0296
0.0300 ± 0.0021% S	0.2511	0.0415	0.0295
	0.2515	0.0420	0.0299
	Avg =	0.0422	0.0297
	s =	0.0006	0.0002

*Results based upon linear, forced through origin calibrations utilizing JK 21 Low Alloyed Steel @ 0.1741% C and EURO 578-1 Ferro-Molybdenum @ 0.065% S.

Procedure / Method 2

1. Prepare the instrument as outlined in the operator's instruction manual.
2. Determine the instrument Blank.
 - a. Login a minimum of three Blank replicates.
 - b. Add ~0.4 g of 501-078 Iron Powder and ~0.6 g of 501-636-HAZ V₂O₅ to a previously heated crucible and thoroughly mix.
 - c. Add ~1.5 g of 763-266 LECOCEL to crucible, covering the accelerators.
 - d. Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
 - e. Repeat steps 2b through 2d a minimum of three times.
 - f. Set the blank by following the procedure outlined in the operator's instruction manual.
3. Calibrate/Drift Correct.
 - a. Login a minimum of three Standard replicates for each calibration/drift reference material to be used.
 - b. Add ~0.4 g of 501-078 Iron Powder and ~0.6 g of 501-636-HAZ V₂O₅ to a previously heated crucible and thoroughly mix. Tare the crucible and accelerators.
 - c. Weigh ~0.25 g of a suitable calibration/drift reference material into the crucible and thoroughly mix.
 - d. Enter the mass and sample identification into the appropriate replicate fields.
 - e. Add ~1.5 g of 763-266 LECOCEL to crucible, covering the reference material and accelerators.
 - f. Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
 - g. Repeat steps 3b and 3f a minimum of three times for each calibration/drift reference material intended for calibration/drift.
 - h. Calibrate/drift correct by following the procedure outlined in the operator's instruction manual.
 - i. Verify the calibration by analyzing ~0.25 g of another suitable reference material following steps 3b through 3f and confirm that the results are within an acceptable tolerance.
4. Sample Analysis.
 - a. Login a sample with an appropriate number of replicates.
 - b. Add ~0.4 g of 501-078 Iron Powder and ~0.6 g of 501-636-HAZ V₂O₅ to a previously heated crucible and thoroughly mix. Tare the crucible and accelerators.
 - c. Weigh ~0.25 g of a ferroalloy sample into the crucible and thoroughly mix.
 - d. Enter the mass and sample identification into the appropriate replicate fields.
 - e. Add ~1.5 g of 763-266 LECOCEL to crucible, covering the sample and accelerators.
 - f. Place the crucible on the furnace pedestal (or appropriate autoloader position if applicable) and initiate analysis.
 - g. Repeat steps 4a through 4f as necessary.

Typical Results / Method 2*

Sample	Mass (g)	% Carbon	% Sulfur
BCS 203/5	0.2514	0.0421	0.0162
Low Carbon	0.2563	0.0427	0.0164
Ferro-Chromium	0.2542	0.0447	0.0164
0.043% C	0.2625	0.0428	0.0160
	0.2567	0.0437	0.0163
	Avg =	0.0432	0.0163
	s =	0.0010	0.0002
Ferro-Chromium	0.2544	0.0337	0.0008
	0.2547	0.0328	0.0008
	0.2524	0.0349	0.0008
	0.2540	0.0341	0.0008
	0.2521	0.0335	0.0008
	Avg =	0.0338	0.0008
	s =	0.0008	< 0.0001
EURO 584-1	0.2512	0.0423	0.0292
Ferro-Titanium	0.2506	0.0432	0.0294
0.0445 ± 0.0043% C	0.2512	0.0426	0.0291
0.0300 ± 0.0021% S	0.2538	0.0437	0.0295
	0.2534	0.0421	0.0301
	Avg =	0.0428	0.0294
	s =	0.0006	0.0004

*Results based upon linear, forced through origin calibrations utilizing JK 21 Low Alloyed Steel @ 0.1741% C and EURO 578-1 Ferro-Molybdenum @ 0.065% S.

