



# Quantitation of Four Regulated Phthalates in Plastic Consumer Products by Gas Chromatography– Time-of-Flight Mass Spectrometry

LECO Corporation; Saint Joseph, Michigan USA

Key Words: Phthalates, GC/MS, GC-MS, Polymer Additives

### 1. Introduction

Phthalates are a class of chemical compounds commonly used to soften and add pliability to plastic materials. Over the past decade, phthalates have come under increased scrutiny as concerns about their potentially harmful effects on human health have arisen. In 2015 the European Union announced a ban on four specific phthalates commonly used in consumer products such as toys, coated fabrics, building materials, and office supplies. The ban included the following phthalates: diisobutyl phthalate, dibutyl phthalate, benzylbutyl phthalate, and bis(2-ethylhexyl) phthalate. This application note demonstrates the use of LECO's Pegasus<sup>®</sup> BT GC-TOFMS system for quantitation of these phthalates in a range of commercial products including two children's toys, an ink pen barrel, and cross-linked high-density polyethylene (CL-HDPE) pipe used in residential plumbing applications.




Figure 1. Total Ion Chromatogram (TIC) and Extracted Ion Chromatogram (XIC m/z 149.05) Showing Detection of 3 out of 4 Regulated Phthalates in Cross-linked High-Density Polyethylene water pipe.

## 2. Experimental

Four different plastic products were tested for the regulated phthalates: two children's toys, an ink pen barrel, and pipe used in residential plumbing. The goal was to detect the phthalates present on the surface of these materials; therefore exhaustive extraction techniques requiring sample grinding prior to extraction were not employed. For each sample, a 1 to 2 gram portion was placed in a 15 mL scintillation vial and the weight recorded. A 10 mL aliquot of chloroform was added, and the vial was placed in an ultrasonic bath for 20 minutes. After sonication, a 1mL portion of each extract was transferred to GC autosampler vials for analysis. Calibration standards for each of the phthalates were prepared in chloroform at a concentration range from 20 to 10,000  $pg/\mu$ L. Plastic extracts and standards were analyzed using the instrument conditions shown below in Table I. External standard calibration curves for each of the phthalates were generated within the ChromaTOF<sup>®</sup> brand software, and were used to quantitate the levels of each regulated phthalate in the extracts.

Table I. GC-TOFMS (Pegasus BT) Conditions	
---	--

Gas Chromatograph	Agilent 7890 with L-PAL3 Autosampler		
Injection	1µL, split 10:1 @ 300 °C		
Carrier Gas	He @ 1.4 mL/min, Constant Flow		
Column	Rxi-5ms, 30 m x 0.25 mm i.d. x 0.25 $\mu$ m coating (Restek, Bellefonte, PA, USA)		
Oven Program	50 °C (1 min), to 330 °C @ 30 °C/min (20 min)		
Transfer Line	300 °C		
Mass Spectrometer	LECO Pegasus BT		
Ion Source Temperature	250 °C		
Mass Range	40-700 m/z		
Acquisition Rate	10 spectra/s		

### 3. Results and Discussion

Table II shows the quantitative results obtained from the surface extraction of the four plastic products. The values are expressed in ng/g.

Sample ID	Diisobutyl phthalate	Dibutyl phthalate	Benzyl butyl phthalate	Bis(2-ethylhexyl) phthalate	
CL HDPE Pipe	491	817	n.d.	1011	
Ink Pen Barrel	2741	5077	973	114450	
Children's Toy #1	3728	5457	404	6894	
Children's Toy #2	1858	1934	619	9561	
Results shown in ng/g					

The GC-MS method used for this analysis demonstrated excellent linearity across the entire calibration range for all four of the phthalates measured. Table III shows the phthalates and their corresponding linearity values. Figure 2 shows a representative calibration plot and spectra for Bis(2-ethylhexyl) phthalate.

Table III. Phthalates and their corresponding linearity values

Phthalate	R <sup>2</sup> Value
Diisobutyl phthalate	0.99999
Dibutyl phthalate	0.99995
Benzyl butyl phthalate	1
Bis(2-ethylhexyl) phthalate	0.99999




Figure 2. GC-MS Calibration Curve and Mass Spectra (Deconvoluted and NIST Library match) for Bis(2-ethylhexyl) Phthalate.

A representative GC-MS total ion chromatogram (TIC) and extracted ion chromatogram (XIC) of m/z 149.05 for the regulated phthalates in Children's Toy #1 are shown in Figure 3.




Figure 3. GC-MS Chromatograms for Toy #1 Extract.

Since the Pegasus BT offers maximum sensitivity and reproducibility with full mass range data acquisition, the same analyses used for the targeted phthalate quantitation can also be interrogated for other untargeted compounds of interest. A representative example is shown in Figure 4 where several additional polymer additive-type compounds were identified.




Figure 4. Example showing the automated identification of several non-targeted compounds within the CL HDPE Extract.

Another important attribute of the system is NonTarget Deconvolution<sup>®</sup> (NTD<sup>®</sup>). Non-Target analysis is a rapidly expanding area of importance in analytical sciences. The *Pegasus* BT was built from the ground up to support these types of investigations. *NTD* is an automated function within the *ChromaTOF* brand software and very effectively identifies known unknowns (known to be in the NIST library but not known to be in the sample a priori). The *Pegasus* BT, being a time-of-flight mass spectrometer, provides data perfectly suited for mass spectral deconvolution algorithms (note scanning instruments are challenged due to spectral skewing). The example shown in Figure 5 demonstrates the ability of *ChromaTOF's* NTD to effectively identify two analytes that were not separated chromatographically. This was done in a fully automated fashion and resulted in excellent spectral similarity scores against the NIST 2017 database. Bis(2-ethylhexyl)Phthalate and 2,4-bis(1-methyl-1-phenylethyl)-Phenol were identified with similarity scores of 901/1000 and 863/1000, respectively.




Figure 5. Example demonstrating effective deconvolution of Bis(2-ethylhexyl)Phthalate and 2,4-bis(1-methyl-1-phenylethyl)-Phenol.

## 4. Conclusion

This application note effectively demonstrates the ability of the *Pegasus* BT GC-TOFMS system to quantitate four regulated phthalates in plastic consumer products. The Time-of-Flight MS system demonstrated excellent sensitivity and linearity for this application. The ability to quantitate targeted analytes while still maintaining the ability to identify other relevant compounds within the same analytical data file shows the *Pegasus* BT GC-TOFMS to be a powerful tool for not only routine screening of regulated components, but also additional characterization of consumer product formulations.





LECO, Pegasus, ChromaTOF, and StayClean are trademarks of LECO Corporation.

LECO Corporation | 3000 Lakeview Avenue | St. Joseph, MI 49085 | Phone: 800-292-6141 | 269-985-5496 info@leco.com • www.leco.com | ISO-9001:2008 | HQ-Q-994 | LECO is a registered trademark of LECO Corporation.