

Determination of 277 Pesticides at Trace Levels by GC–TOFMS

Sjaak de Koning, Monika Gumpendobler • LECO Instrumente, Mönchengladbach, Germany

Key Words: TOFMS, Pesticides, Trace Analysis, Quantification, Food

1. Introduction

The determination of trace contaminants in complex matrixes, such as food, often requires extensive sample extraction and preparation regimes prior to instrumental analysis. In recent years, the use of mass spectrometry in instrumental analysis has enabled the 'extraction' of an analyte at the detection stage of a method. This has led to the development of crude extractions—due to the need for rapid and universal methods, covering a large number of compounds down to the levels required by legislation [1,2]. For example, the QuEChERS method for pesticides analysis [3], which is based on dispersive solid phase extraction (DSPE), uses lower volumes of solvent and little or no clean-up.

The purpose of this study is the evaluation of a GC-TOFMS method for the trace analysis of pesticides in food stuffs.

2. Instruments and Methods

In this study, measurements were made with a LECO Pegasus® GC–TOFMS system. This system consists of an Agilent 6890N gas chromatograph with a LECO Pegasus Time-of-Flight Mass Spectrometer (TOFMS) as a detector. For all analyses in this study, the analytical column was a 30 m x 0.25 mm ID x 0.25 μ m df DB-5ms. Helium was used as a carrier gas at a constant flow of 1.4 mL/minute. The transfer line to the TOFMS consisted of the last 20 cm of the analytical column. The detector was optimized to have a signal-to-noise ratio >1000 for a 200 pg/ μ l injection of hexachlorobenzene. Data are acquired with a mass range of 50 to 550 amu, an acquisition rate of 10 spectra/second at a detector voltage of 1850 V. In this one-dimensional separation, the temperature program of the GC was as follows.

Rate (°C/min)	Target Temp (°C)	Duration (min)
Initial	95	1.50
20	190	0.00
5	230	0.00
25	300	20.00

For all quantitative work, pesticides standards in the range of 2.5 to 1000 pg/ μ l were prepared in acetonitril. Injections of 2 μ l were performed by a CTC CombiPAL with a Gerstel CIS4 PTV injector in splitless mode at a temperature of 40°C. After 0.25 minutes the injector was heated to 280°C (held for 8 minutes) at a rate of 12°C/second.

3. Results

First, after analysis of the analytical standards (n=7), for all pesticides the average retention time was determined for 500 pg/ μ l (see Figure 1 for TIC chromatogram). This standard was used as a reference standard for further work. For all pesticides the analytical data as determined are shown in Table 1.

As can be expected, many coelutions do occur. Figure 2 shows an example of how the True Signal Deconvolution™ algorithm of the ChromaTOF® software can be used to identify coeluting peaks, delivering true peak spectra.




Figure 1. TIC chromatogram showing all pesticides studied.




Figure 2. Zoomed TIC and EIC of coeluting compounds. Peak identifications are displayed at the peak markers.

4. Conclusions

It is shown in this work that the LECO Pegasus GC–TOFMS equipped with a Gerstel CIS4 PTV injector is able to analyze hundreds of pesticides at trace levels (down to 10 ppb) in a single run. This work is advantageous to the food and flavor industry where pesticide contamination poses severe health risks. The ability to acquire full mass range spectra allows ChromaTOF software to perform True Signal Deconvolution, enabling peak finding, identification, and quantification of overlapping pesticides.

5. References

- [1] Commission Directive 2003/13/EC, *Off. J. Eur. Com.* L41 (2003) 33.
- [2] Commission Directive 2003/14/EC, *Off. J. Eur. Com.* L41 (2003) 37.
- [3] M. Anastassiades, S.J. Lehotay, D. Stajnbaher, F.J. Schenk, *J. AOAC Int.* 86 (2003) 412.

Table 1. Pesticides analytical data.

Analyte	Compound Name	LOD* (pg/µl)	Absolute R.T. (s)	Quant Mass	Analyte	Compound Name	LOD* (pg/µl)	Absolute R.T. (s)	Quant Mass
1	o-Cresol	1.07	184.594	108	81	Triethazine	5.02	491.479	200
2	p-Cresol	5.21	192.779	107	82	Terbutylazine	1.79	493.394	214
3	m-Cresol	2.29	193.294	107	83	Cyanophos	3.84	493.979	109
4	Omethoate	21.49	203.979	110	84	Propyzamide	5.88	496.879	173
5	Aniline, 2,5-dimethyl-	17.94	233.549	106	85	Diazinon	1.69	498.479	137
6	2,5'-Dimethylaniline	0.71	234.294	121	86	Fonofos	1.05	498.879	109
7	m/p-Chloroaniline	2.10	250.479	127	87	Pyrimethanil	1.01	505.394	198
8	Methamidophos	25.56	258.879	94	88	Dimethirimol	13.58	505.858	166
9	Dichlorvos	32.72	264.779	109	89	Disulfoton	14.40	512.294	153
10	Dichlorfos	16.67	264.894	185	90	Terbacil	3.37	513.694	161
11	Bendiocarb	3.38	276.779	126	91	Fenfuram	5.68	515.994	201
12	Dichlobenil	1.70	308.979	171	92	Triallate	19.36	520.379	268
13	ETPC	1.74	312.994	128	93	d-HCH	14.92	521.394	181
14	3,5-Dichloroaniline	4.64	334.479	161	94	Pirimicarb	1.86	526.094	166
15	Mevinphos-trans	9.89	339.979	127	95	Tridemorph	43.30	532.079	128
16	Ethiofencarb - M	3.44	342.094	168	96	Formothion	106.89	535.479	125
17	3,4-Dichloraniline	6.68	345.979	161	97	Metobromuron	23.54	538.079	170
18	Chlormephos	2.31	346.879	121	98	Pentachloroaniline	9.41	539.479	265
19	Vernolate	1.27	349.194	128	99	Endosulfan ether	1.24	541.879	69
20	Etridiazole	1.09	351.694	211	100	Phosphamidon (Z)	2.45	542.694	127
21	Propham	3.11	354.194	179	101	Dichlofenthion	9.50	546.279	223
22	Bromophos-ethyl	35.66	356.394	242	102	Dimethachlor	2.25	546.994	197
23	Methacrifos	5.77	357.394	180	103	Chlorpyrifos-methyl	3.60	553.994	286
24	α-HCH	1.50	357.694	181	104	Metribuzin	6.36	554.179	198
25	Methacrifos:2	1.53	368.294	180	105	Vinclazoline	9.65	558.179	124
26	Crimidine	2.34	374.694	171	106	Methyl parathion	22.01	562.094	263
27	Tebuthiuron	14.19	379.479	156	107	Alachlor	12.84	562.794	188
28	o-Phenylphenol	2.90	379.979	170	108	Tolclofos-methyl	4.05	563.379	265
29	Benzene, pentachloro-	0.55	380.394	250	109	Fuberidazole	22.74	567.494	184
30	β-HCH	1.10	385.894	181	110	Carbaryl	48.06	570.079	144
31	Methiocarb	1.29	389.594	168	111	Metalaxyl	8.65	571.479	132
32	TEPP	149.99	391.979	161	112	Ametryn	16.48	572.394	227
33	Phenol, 2,3,4,5-tetrachloro-	25.45	392.594	232	113	Heptachlor	9.20	574.094	272
34	ETU	75.27	393.194	102	114	Fenchlorphos	4.18	575.394	285
35	Methomyl	56.87	394.479	105	115	Paraoxon	17.96	577.179	109
36	Heptenophos	7.67	397.279	124	116	Demeton-S-methylsulfone	466.25	584.179	169
37	Chlorfenprop-methyl	4.06	407.779	165	117	Propargit-Metabolit	5.03	585.579	135
38	Tecnazene	2.70	408.294	261	118	Pirimiphos methyl	15.26	589.479	180
39	Thionazin	2.07	411.294	143	119	Terbutryn	5.41	592.579	170
40	Propoxur	217.64	412.894	110	120	Fenitrothion	18.13	593.794	277
41	Propachlor	0.99	413.894	120	121	Ethofumesate	6.61	596.279	161
42	Demeton-O	20.13	413.979	115	122	Bromacil	0.62	598.179	205
43	Demeton-S-methyl	20.19	419.779	109	123	Linuron	8.80	601.679	61
44	Fluometuron	4.96	420.879	72	124	Malathion	29.16	605.494	173
45	Diphenylamine	0.84	423.294	169	125	Metolachlor	9.90	611.979	238
46	Tributyl phosphate	1.21	425.094	155	126	Chlorpyrifos-ethyl	12.07	614.894	199
47	Ethoprosfos	3.41	425.179	158	127	Diethofencarb	12.26	617.379	168
48	Cycloate	1.62	427.194	154	128	Benthiocarb	1.69	617.894	100
49	2,4,5-Trimethacarb	72.69	433.179	121	129	Cyanazine	37.01	620.794	225
50	Trifluralin	1.38	434.594	264	130	Fenthion	3.67	621.394	278
51	Chlordimeform	6.77	435.579	181	131	Anthraquinone	7.99	624.394	180
52	Dicrotophos	4.39	436.194	127	132	Parathion-ethyl	23.42	626.294	291
53	Dichlofluanid-Metabolit (DMSA)	25.86	436.379	200	133	Triadimefon	15.81	630.279	208
54	Sulfotep	7.30	440.579	322	134	Carbetamide	138.82	632.979	119
55	Monocrotophos cis/trans	21.65	443.179	127	135	Chlorthion	15.50	636.679	109
56	Diallate cis	2.97	451.794	234	136	Nitrothal-isopropyl	7.95	636.794	212
57	Phorate	2.69	452.494	121	137	Telodrin	1.88	638.894	103
58	Diallate trans	4.38	460.594	234	138	Trichloronat	5.03	639.294	269
59	Hexachlorobenzene (HCB)	0.87	461.994	284	139	Crufomate	28.43	639.679	182
60	Thiometon	4.54	464.179	125	140	Butralin	9.31	639.994	266
61	Anisole, 2,3,4,5,6-pentachloro-	5.67	466.779	237	141	Fenson	3.94	645.279	141
62	Demeton-S	2.64	467.794	88	142	Pirimiphos ethyl	12.63	645.794	168
63	Dimethoate	8.26	469.094	125	143	Dodemorph-cis	7.70	649.294	154
64	Dichloran	4.23	469.794	206	144	Isodrin	4.19	662.479	195
65	2,3,5-Trimethacarb	58.46	469.879	121	145	Metazachlor	2.44	665.694	132
66	Ethoxyquin	10.86	471.694	202	146	Penconazole	7.60	671.879	159
67	Carbofuran	340.51	473.079	164	147	Endosulfan sulfate	13.92	673.779	239
68	Simazine	18.86	474.979	201	148	Chlzolinate	32.10	674.394	188
69	Atrazine	7.30	479.179	200	149	Isofenphos	7.99	676.279	213
70	Chlorbufam	31.05	479.594	127	150	Oxychlordan	7.12	676.694	154
71	Monolinuron	7.40	480.594	126	151	Heptachlor epoxide	22.37	677.094	353
72	Dimethipin	78.15	481.179	118	152	trans-Chlorfenvinphos	17.07	678.594	267
73	Propazine	3.27	482.694	172	153	Mecarbam	29.30	682.194	131
74	Cyromazine	42.38	483.494	151	154	Bioallethrin	16.58	684.794	123
75	Quintozen	14.51	485.479	142	155	Phenthoate	8.59	686.779	125
76	Benazolin	31.05	487.594	170	156	Pyrethrin	15.35	687.079	123
77	Tolyfluanid-Metabolit (DMST)	37.88	488.494	214	157	Quinalphos	5.49	688.279	146
78	Dioxathion	85.19	488.894	270	158	Tiabendazole	42.16	689.394	174
79	Propetamphos	2.03	490.894	194	159	Procymidone	14.09	692.479	283
80	g-HCH	0.81	491.194	181	160	Triflumizole	6.84	696.179	73
					161	Crotoxyphos	76.27	698.779	166
					162	Triadimenol	11.68	705.679	112

Analyte	Compound Name	LOD* (pg/ μ l)	Absolute R.T. (s)	Quant Mass	Analyte	Compound Name	LOD* (pg/ μ l)	Absolute R.T. (s)	Quant Mass
163	Bromophos-ethyl	7.00	709.794	145	228	Propargite II	56.37	923.979	201
164	Chlordane , trans	3.96	713.394	125	229	Diflufenicam	2.45	924.294	266
165	Chinomethionat	40.97	715.379	206	230	Piperonyl butoxide	4.38	930.379	176
166	o,p'-DDE	7.19	717.179	246	231	Benzoylprop ethyl	1.16	939.594	105
167	Tetrachlorvinphos	30.55	719.679	329	232	Pyridafenthion	47.79	945.879	340
168	Vamidothion	60.67	721.094	145	233	Tetramethrin	21.15	951.194	164
169	Paclbutrazide	16.68	722.579	236	234	Phosmet	1.83	951.994	160
170	Endosulfan I	14.02	733.579	195	235	Bifenthrin	1.03	956.794	181
171	Ditalimfos	6.12	736.779	130	236	Bromopropylate	10.35	956.994	341
172	Diethatyl ethyl	10.62	737.479	188	237	Tetramethrin:2	2.64	959.494	164
173	Flutriafol	3.24	741.179	123	238	Methoxychlor	3.26	961.979	227
174	Bromfenvinphos	16.83	741.879	267	239	Fenpropathrin	33.43	965.279	181
175	Fenamiphos	13.61	744.779	154	240	Bifenox	132.70	970.879	314
176	Jodfenphos	10.68	750.594	377	241	Furathiocarb	353.79	977.379	194
177	Chlorfenson	12.91	752.894	175	242	Phenothrin	2.85	977.794	183
178	Imazalil	14.26	753.479	215	243	Tetradifon	10.23	979.279	159
179	Prothiofos	26.34	756.894	309	244	Phenothrin:2	1.69	982.994	183
180	Profenofos	92.36	764.594	208	245	Phosalone	2.23	984.694	182
181	DDE p,p'	9.10	770.894	246	246	Leptophos	1.61	985.294	171
182	Flamprop-methyl	1.55	775.579	105	247	Azinphos-Methyl	14.15	987.694	160
183	Myclobutanil	5.82	777.879	179	248	Mirex	1.38	1001.89	272
184	Dieldrin	11.65	777.894	263	249	Fenarimol	11.85	1006.98	221
185	o,p'-DDD	7.02	781.279	235	250	Azinphos-ethyl	5.52	1012.48	160
186	Bupirimate	19.48	783.579	273	251	Dialifos	17.69	1014.98	208
187	Buprofezin	9.97	784.094	172	252	Bitertanol B	1.02	1028.39	170
188	Diclobutrazol	6.74	787.079	147	253	Permethrin	1.09	1031.19	183
189	Cinerin	40.88	790.779	123	254	Bitertanol B:2	24.32	1032.59	248
190	Binapacryl	366.26	798.379	83	255	Coumaphos	13.12	1035.29	226
191	Nitrofen	28.46	812.294	283	256	Permethrin:2	3.66	1035.99	183
192	Endrin	40.76	812.794	263	257	Pyridaben	2.62	1037.18	147
193	Fluazifop-butyl	10.97	818.194	282	258	Prochloraz	11.09	1037.69	180
194	Chlorobenzilate	2.38	828.794	139	259	Cyfluthrin	32.50	1051.68	163
195	Fenthion sulfoxide	40.22	830.994	279	260	Cyfluthrin:2	33.05	1056.18	163
196	Flamprop-methyl:2	1.40	831.279	105	261	Cyfluthrin:3	164.10	1058.38	163
197	Fensulfothion	18.58	831.994	293	262	Cyfluthrin:4	126.99	1059.98	163
198	Endosulfan II	22.39	832.594	195	263	Cypermethrin	95.83	1065.28	181
199	Etaconazole cis	16.93	833.779	173	264	Cypermethrin:2	123.77	1069.78	181
200	Oxadixyl	5.57	839.979	132	265	Cypermethrin:3	106.83	1071.98	181
201	Etaconazole trans	11.95	840.379	173	266	Quinofop-ethyl	13.51	1073.08	299
202	DDD p,p'	3.83	843.894	235	267	Fluvalinate	249.14	1076.68	250
203	Ethion	7.12	845.579	231	268	Decachlorobiphenyl (PCB 209)	17.30	1084.28	498
204	DDT p,p'	1.23	846.994	235	269	Fenvalerate à	147.63	1106.08	167
205	Chlorthiophos iii	11.57	849.479	269	270	Fenvalerate à	463.74	1106.39	225
206	Chlorthiophos	9.04	849.594	269	271	Fluvalinate:2	209.51	1114.38	250
207	Benodanil	10.92	859.894	231	272	Esfenvalerate	207.61	1116.38	167
208	Tetrasul	14.03	867.579	252	273	Fenvalerate à:2	448.34	1116.64	225
209	Triazophos	12.50	872.094	161	274	Difenoconazole cis	14.59	1129.89	323
210	Sulprofos	5.16	872.679	156	275	Difenoconazole trans	18.47	1133.19	323
211	Pyrethrin	55.01	878.779	123	276	Deltamethrin	74.69	1135.23	181
212	Benalaxyl	2.26	882.894	148	277	Deltamethrin:2	260.29	1146.22	181
213	Carbophenothion	4.10	887.379	157					
214	Norflurazon	12.97	889.079	145					
215	Cyanophenphos	5.32	889.879	157					
216	Edifenphos	19.54	890.379	109					
217	Propiconazole cis	15.78	892.479	173					
218	Propiconazole cis	50.68	892.694	259					
219	Lenacil	13.62	896.794	153					
220	Chloridazon	45.17	899.394	221					
221	Propiconazole trans	5.48	900.379	173					
222	Propiconazole trans	42.06	900.594	259					
223	o,p'-DDT	2.71	901.694	235					
224	Hexazinone	2.69	907.594	171					
225	Nuarimol	4.39	915.079	139					
226	Tebuconazole	16.70	918.179	250					
227	Propargite I	34.49	922.073	201					

*LOD is expressed as three times the noise.



LECO Corporation • 3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977
 info@leco.com • www.leco.com • ISO-9001:2000 • No. FM 24045 • LECO is a registered trademark of LECO Corporation.