



Brinell Indentation Microscope

The KingScan IV Brinell Measuring Microscope System is completely automated, providing consistent and accurate results down to .01mm and is ASTM E-10 and ISO compliant. In a fraction of a second, it measures indentation diameters and instantaneously transfers the results to your computer monitor.

New-digital camera is faster than previous models

Works with any off the shelf USB 2.0 cable

Direct computer connection — no junction box needed

Windows 7, Vista and XP, 32- and 64-bit compatible

MINIMUM REQUIREMENTS:

1.8 Gig Celeron CPU or Higher
256 MB DDR Ram or More
Intel 3D onboard Graphics /
Graphics N-VIDIA G-Force 3 or Higher
1 - USB 2.0 Port
Windows XP Service Pack 1 or 1a,
Windows 2000 Service Pack 3 (or Higher)

RECOMMENDED

2.4 Gig P4 CPU or Higher
512 MB 333 - 400 MHZ Dual Channel
DDR Ram or Higher
On-board Intel PCI Express Graphics /
N-VIDIA G-Force 5 or Higher
1 - USB 2.0 Port
Windows 7 or Windows Vista