

# LCB3100 Brinell

## Load Cell Hardness Tester

Using state-of-the-art load cell technology, this is one of the most versatile and easy-to-use Brinell Hardness Testers available. Operation is handled through a user-friendly touch-pad, virtually eliminating errors by removing the use of knobs and the stacking of weights typically associated with most testers. Simply choose the load on the operation panel and let the tester do the rest.

- Twelve selectable loads to choose from
- Heavy-duty cast housing with enclosed load cell
- Direct load method with load-cell feedback
- Auto-apply, -hold, and -removal of loads
- Selectable dwell times
- User-friendly touch-pad operation panel

