

# MX400

## Mounting Press

Our fully automatic mounting press, the MX400, is your solution for faster, more sophisticated mounting. The MX400 boasts a slim, 8-inch wide profile to maximize lab space. A touch-screen display provides easy, user-friendly operation.

### Features

- Novel solid steel crossbar closure design is durable, easy to clean, avoids jamming, and allows for quick, one-handed operation
- Intuitive swipe-and-touch user interface makes it easy for you to select heating, cooling, and pressure parameters
- Create fully customizable methods or use a LECO-developed method to fit your exact needs
- Variable modes for mounting delicate samples using all common types of mounting materials
- Low thermal mass for rapid heating and cooling rates for superior productivity
- Reduce mounting bottlenecks and increase throughput with full size dual mount capability in 1.25 inch, 1.50 inch, 30 mm, 40 mm, and 50 mm mold size



**Easily mount samples  
in four simple steps.**

**LECO**  
EMPOWERING RESULTS