
## Analysis of Pesticides in Edible Oils by GC×GC–TOFMS

European Technical Center
LECO Instrumente GmbH, Mönchengladbach, Germany



Peak markers indicate a wide range of organochlorine and organophosphorus pesticides that were detected with an s/n ranging from 100 to 1000.

## **Columns:**

30 m x 0.25 mm x 0.25  $\mu$ m VF-5ms 2 m x 0.1 mm x 0.2  $\mu$ m VF-17ms

## TOFMS:

50 to 550 m/z @ 100 spectra/s

For questions on this analysis e-mail us at: life\_science@leco.com

For a complete listing of snapshots and application notes visit us on the web at www.leco.com



